

# VIVASIGHTDL

By ETVIEW



Left-sided VIVASIGHT-DL™

# VIVA SIGHT DL

THE ONLY PLATFORM THAT ALLOWS FOR CONTINUOUS, REAL-TIME, MANAGEMENT OF THE AIRWAY,  
NOW AVAILABLE IN SIZES: 35F, 37F, 39F AND 41F.

VivaSight-DL's accurate positioning and continuous control overcome the limitations of current procedures with a bronchoscope

VivaSight-DL is the world's first sterile, single-use, double-lumen airway tube with an integrated high-resolution camera. VivaSight-DL overcomes the placement and positioning challenges of current DLTs, and offers unparalleled continuous visualization throughout the procedure.

| Code    | Description        | Effective ID (mm) | Bronchial OD (mm) |
|---------|--------------------|-------------------|-------------------|
| DLVT35L | 35 Fr VivaSight-DL | 4                 | 10.5              |
| DLVT37L | 37 Fr VivaSight-DL | 4.2               | 11.5              |
| DLVT39L | 39 Fr VivaSight-DL | 4.8               | 11.5              |
| DLVT41L | 41 Fr VivaSight-DL | 5                 | 12                |

ETVIEW MEDICAL Ltd. | The Trendlines Building  
Misgav Business Park | M.P. Misgav 20174 | Israel  
T. +1 567 248 4822 | [www.etview.com](http://www.etview.com)

**Caution:** Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner.  
DL clearance/certification: FDA, CE (MDD 93/42/EEC), ISO 9001:2000, ISO 13485:2003  
VivaSight-DL is a registered trademark of ETView Medical Ltd.